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 Joinstar Biomedical Technology Co.,Ltd.
Address:10th Floor ,Administration Building,
NO.519,XingGuo RD.,Yuhang Economic and
Technological Development Zone, Hangzhou ,
Zhejiang, China, 311188
E-mail: market@joinstar.cn
Tel: 0086-571-89023160

 Lotus NL B.V.
Address: Koningin Julianaplein 10, 1c Verd, 2595AA,
The Hague, Netherlands.
E-mail: peter@lotusnl.com
Tel: +31644168999

COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device

Instructions For Use

【PRODUCT NAME】

COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device

【SPECIFICATIONS】

1 test/kit, 5tests/kit, 10tests/kit, 25tests/kit.

【INTENDED USE】

This kit is used to detect the presence of SARS-CoV-2 antigens in human anterior nasal swab specimens. It can be used for self-monitoring in case of symptoms characteristic of COVID-19, however, it does not replace tests based on nucleic acid detection (for example PCR analysis).

A positive result indicates that the sample contains SARS-CoV-2 antigens. A negative result does not rule out the possibility of infection.

This kit is intended for out-of-laboratory use, self-monitoring. The test results from the kit are for information purposes only and are not a substitute for tests performed by healthcare professionals. An accurate diagnosis can only be made based on the patient's clinical symptoms and the results of other laboratory tests.

The complete test (including all components) can be used only once, after use place the devices in the waste container as indicated below.

Note: This is a self-test. The results obtained are for information only. If you experience COVID-19 symptoms, avoid personal contact and follow the regulations in force.

【TEST PRINCIPLE】

The COVID-19 (Colloidal Gold) antigen assay uses the principle of lateral flow immunoassay to detect the presence of SARS-CoV-2 virus nucleocapsid protein in nasal swab specimens. When the sample is dropped into the sample well, if the SARS-CoV-2 virus nucleocapsid protein is present, it is bound to color-labeled antibodies that can uniquely recognize the protein. The sample then enters the test area of the membrane due to the capillary action. There are two bands in this region, the first (labeled T) contains antibodies that can uniquely

recognize the SARS-CoV-2 virus nucleocapsid, and the second band (labeled C) binds only antibodies that recognize color-labeled antibodies to the membrane. Thus, the T-labeled band may bind to SARS-CoV-2 virus nucleocapsid proteins that have previously bound a colored antibody, so the T-label only shows a color reaction if this viral protein was present in the nasal swab sample. If no viral protein was present in the sample, the T mark would not discolor the band. In the C-band, the color-labeled antibodies bind, so a color reaction must take place in each case. If the C-labeled band does not stain, the color-labeled antibodies did not reach the test surface. If there is no color reaction here, the test cannot be evaluated.

【COMPONENTS】

1. Test cassette
2. Extraction buffer
3. Sterile swab
4. Instructions for use

A timer is needed, which is not included in the kit.

Note: Components from different batches may not be mixed.

Mixing different batches of equipment will invalidate the test. Unpack the devices only at the step indicated below.

【STORAGE AND SHELF LIFE】

The test is stable for 2 years if all parts are stored in sealed pouches at +2°C to 30°C and the packaging of all devices is unopened and undamaged.

The test should be stored in the sealed pouch until use. Do not freeze the test and do not use it after the expiry date.

Please refer to the package for the date of manufacture next to the mark  and expiry date next to the mark .

【SAMPLE REQUIREMENTS】

Read the entire leaflet before performing the test.

Collection of samples:

1. Please disinfect your hands thoroughly before use!
2. Then open the package containing the sampling wand and remove the swab from the packaging.
3. Gently insert the entire soft end of the swab about 2 cm into one nostril.


4. Using medium pressure, slowly rub the swab in a circular motion around the inside of your nostril 4 times for a total of 15 seconds.
5. Repeat the same procedure with the same swab in the other nostril.

Washing the sample into the extraction liquid:

1. Open the tube containing the extraction buffer.
2. The soft end of the sampling swab is dipped below the liquid level of the sample extraction buffer after sampling.
3. In the solution, rotate the swab and squeeze the side of the tube, repeating this three times. The soft end of the sampling swab should be kept in the solution for at least 15 seconds.
4. Then press the soft end of the swab three times against the inner wall of the tube.
5. Then place the swab in the waste bag.
6. Close the tube containing the buffer solution with the dropper cap.


Storage of samples:

Drop the sample onto the test cassette as soon as possible, but no later than one hour after the previous process. If you do not perform the test for more than an hour after sampling, the test result will be considered invalid.

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1. The test should be performed at room temperature (15 ° C-30 ° C). Before performing the test, make sure that both the test cartridge and the extraction buffer are at room temperature.
2. Remove the test cassette from the sealed pouch and place the test cassette on a flat surface.
3. Place 2 drops of the sample-containing liquid prepared as described in the previous step into the sample hole (labeled S).

4. Please set a timer to 15 minutes. The result can be read 15 minutes after dropping the sample-containing liquid on the cassette.

The check line marked C may appear earlier. However, be sure to wait 15 minutes before reading. The results you read before this are invalid. If no result is available after 20 minutes, it is recommended to repeat the test with a new test cassette.

【SAMPLE PROCESSING AND CLEANING】

1. Place the test cassette, sample extraction buffer and used swab in a waste bag and seal the bag.

2. Throw the bag in the bin.

3. Re-apply the hand sanitizer.


【INTERPRETATION OF RESULTS】

Positive (+): Two lines appear (the lines discolor at C and T). One line should always appear in the control line area (C) and another line, in a different colour, should appear in the test line area (T) (see Result 1). NOTE: The intensity of the colour in the test line area may vary depending on the concentration of SARS-CoV-2 nucleocapsid protein in the sample. Therefore, any colour shading in the test line area (T) should be considered positive.

So regardless of the intensity of the discoloration, if the colored line also appears at the C and T markings, the test is positive! If the C mark does not discolor the result is invalid, the test must be repeated.

Negative (-): A colored line appears at the control area (C). No colored line appears in the test line area (T) (see Result 2). **If the C mark does not discolor the result is invalid, the test must be repeated!.**

Invalid: The colored line does not appear at the control area (C) (see result 3). Insufficient sample volume or an inappropriate procedure are the most likely causes of error. Check the instructions for use and repeat the test with a new kit. If the problem persists, stop using the test kit immediately and contact your Hungarian distributor.

A negative result only gives information about your current condition and does not mean that you do not carry the SARS-Cov-2 virus. The result obtained is for information only, SARS-CoV-2 infection can only be completely ruled out by a nucleic acid detection test (eg PCR test).

My result is positive. What should I do? There is a suspicion of infection with Covid-19. Contact a doctor or your local health authority immediately according to current regulations. Follow local guidelines for self-isolation, avoid contact with other people! A confirmatory PCR test must be performed. You can get more information on the website <https://www.who.int/covid-19>.

My result is negative. What should I do? Continue to follow all applicable rules for contact with others and protective measures. Even if the test is negative, an infection may still be present. In all cases of suspicion, repeat the test after 1-2 days, as the coronavirus cannot be reliably detected at all stages of infection.

My result is invalid. What should I do? Repeat the test, consult a doctor or the Covid-19 test center, if the test result is still invalid, notify the Hungarian distributor of the test.

[RESTRICTIONS]

- The results of this test are for information only. The infection should be confirmed by a healthcare professional by comparing this test with other laboratory results, clinical symptoms, epidemiological data, and additional clinical data.
- In the early stages of infection, low virus concentrations may result in a negative test result.
- Test results depend on the quality of sample collection, processing, transport and storage. Errors can lead to false negative results. If contamination is not avoided during sample processing, a false positive result may be obtained.

[SENSITIVITY AND SPECIFICITY]

Method	RT-PCR		Total Results
COVID-19 Antigen Rapid Test	Results	Positive	Negative
	Positive	98	4
	Negative	4	496
Total Result		102	500

Diagnostic sensitivity of the COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device: 96.1% (90.26%-98.92%)*
Diagnostic specificity of the COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device: 99.2% (97.86%-99.92%)*

Overall agreement of the COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device: 98.7% (97.40%-99.42%)*
* 95% confidence interval.

[CROSS-REACTIONS]

The test cassettes were tested for possible cross-reactions with the following pathogens. No cross-reactivity was observed.

HCoV-229E	HCoV-OC43	HCoV-NL63	MERS-CoV
HCoV-HKU1	Human RSV	Human enterovirus	Human rhinovirus
Human metapneumovirus	Mycoplasma pneumoniae	Parainfluenza	Adenovirus
Influenza B virus (Victoria line)	H1N1 (2009) influenza virus	Influenza A virus H3N2	Avian influenza virus H7N9
Influenza B virus (Yamagata series)	Neisseria meningitidis	Seasonal influenza A H1N1	Streptococcus pneumoniae
Staphylococcus aureus			

[INTERFERING SUBSTANCES]

The test cassettes were tested for possible cross-reactions with the following chemicals and drugs. No cross-reactivity was observed.

Interfering substances	Concentration
Aspirin	30ug/dL
Ascorbic Acid	20mg/dL
Ibuprofen	200ug/dL
Bilirubin	60mg/dL
Chloramphenicol	3ug/dL

! [PRECAUTIONS]

- This kit is intended for in vitro diagnostic use only. Please read these instructions carefully before testing.
- Use only the sample swab and extraction buffer provided in the package and do not substitute components from other kits for the specimen extract provided in this kit.
- Strictly follow the instructions in this manual when performing the test.
- Positive and negative predictive values are highly dependent on prevalence. When disease prevalence is low and SARS-CoV-2 is not or only minimally active, positive test results are more likely to represent a false positive result. When disease prevalence is high, false negative results are more likely.
- Compared to a SARS-CoV-2 RT-PCR test, this test is less sensitive when used to detect patient samples within the first five days of symptom onset.
- Waste and excess samples generated during the test may be disposed of with household waste.
- The contents of the desiccant pad included in the test cartridge package are toxic, do not consume.

[FREQUENTLY ASKED QUESTIONS]

- When can I be tested?
You can be tested at any time, whether you have symptoms or not. Please note that the test result is a snapshot and is only valid for that particular moment. Tests must therefore be repeated as required by the relevant authorities.
- What do I need to do to get the most accurate result?
Always follow the manual instructions exactly. Perform the test immediately after taking the sample. Place the drops from the extraction tube only into the hole of the test cassette, two drops at a time. Too many or too few drops may result in a false or invalid test result.
- The test strip is very discoloured. What is the reason for this or what am I doing wrong?
The reason for a clearly visible discolouration of the test strip is that too many drops from the tube containing the extraction buffer have been dispensed from the extraction tube into the sample hole of the test cassette. The test strip can only absorb a limited amount of liquid. If the control line does not appear or the test strip is heavily discoloured, repeat the test with a new test kit following the instructions for use.
- What should I do if I do not see a control line after the test?
In this case, the test result should be considered invalid. Please repeat the test with a new test kit following the instructions for use.
- I don't know how to interpret the results. What should I do?
If you are not sure about the test result, contact your nearest health authority following the instructions of your local authority..
- Can this test cassette be reused or used by more than one person?
No. The entire test (with all its components) is intended for single use only.
- Do I have to follow spacing and hygiene rules if I have a negative result?
Yes. Even in the event of a negative result, the rules of distance and hygiene must be observed..

[EFFECTIVE DATE AND VERSION]

Effective Date: 2021.10.26

Version: 0

Note: Please refer to the table below to identify various symbols.

	Read the instruction manual
	Shelf life until
	Lot
	Catalogue number
	Attention
	Manufacturer
	Date of Production
	Representative in the European Union
	In vitro diagnostics
	Temperature limit
	Do not reuse
	This product fulfils the requirements of the Directive 98/79/EC on in vitro diagnostic medical device
	Tests per kit
	Self-test

Joinstar Biomedical Technology Co.,Ltd.
Address: 10th Floor, Administration Building,
NO.519,XingGuo RD.,Yuhang Economic and
Technological Development Zone, Hangzhou,
Zhejiang, China, 311188
E-mail: market@joinstar.cn
Tel: 0086-571-89023160

Lotus NL B.V.
Address: Koningin Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands.
E-mail: peter@lotusnl.com
Tel: +31644168999